

## Preprocessing Tutorial

This tutorial demonstrates how to preprocess raw BIDS-formatted MEG data using *Brainstorm* to extract necessary time series on which single and multi-channel burst detection can be applied.

The sample data used in this work can be downloaded from [OpenNeuro](#). This BIDS-formatted dataset consists of resting state MEG recordings from 5 healthy participants who are part of the [OMEGA database](#). In this work, we show results for the first of these subjects (sub-0002) but you may wish to apply the tutorials to all subjects (or your own data) to gain a sense of the inter-subject variability.

For convenience we recommend that all tutorial scripts, sample data, and data derivatives be saved in the same local working directory.

### 1. Download and Prepare MEG Data




To begin, follow the steps in the [Brainstorm OMEGA Tutorial](#) up to and including step 6 (Refine the MEG-MRI Registration). For this work, it is only necessary to complete the steps for sub-0002.

At the end of these steps, all data should be loaded into Brainstorm, co-registered, and ready for preprocessing.

### 2. Preprocess MEG Data


Apply bandpass filter (1-50 Hz) and resampling (250 Hz):

- Drag and drop all the subjects (including sub-emptyroom) into the Process1 box and click on [Run]
- Select process: **Import > Import recordings > Convert to continuous (CTF)**: Continuous
- Add process: **Pre-process > Band-pass filter**: High-pass filter at 1Hz, low-pass filter at 50 Hz, 60dB, Process entire file.
- Add process: **Pre-process > resample**: 250 Hz



Artifact cleaning – remove heartbeats and eyeblinks using automate SSP procedure

- In Process1, select all the subject rest recordings (all the subjects, excluding sub-emptyroom)
- Select process: **Events > Detect heartbeats**: channel name ECG, all file, cardiac
- Add process: **Events > Detect eye blinks**: channel name VEOG, all file, blink
- Add process: **Events > Remove simultaneous**: remove cardiac close to blink, delay:250 ms
- Add process: **Artifacts > SSP > Heartbeats**: cardiac, MEG, use existing SSP
- Add process: **Artifacts > SSP > Eye blinks**: blink, MEG, use existing SSP



### 3. Source Estimation

Complete section 9 (Source Estimation) in the [Brainstorm OMEGA Tutorial](#).


The completion of these steps will result in filtered, source-projected MEG data.

### 4. Single-Channel Methods: Extract Single Schaefer Time Series

The single channel burst detection methods are performed on a single time series representing activity in the right motor cortical region. In order to extract such a time series, we must parcellate the data using an anatomical atlas and export the time series at the region of interest as follows:


In Anatomy tab:

- Right click sub-0003 > Add MNI parcellation > Schaefer2018 > Schaefer2018\_200\_17net




In the Functional data tab:

- Double click on the source file to open a source plot
- In the right panel, select the Scout tab then select Atlas > From subject anatomy > Schaefer2018\_200\_17net (MNI-linear)




- This will load the subject-specific Schaefer parcellation
- Close all figures
- Drag sub-0002 to Process 1 box and select Process sources > [Run]
- Select process: **Extract > Scout time series**: all file, Select scouts from volume: Schaefer2018\_200\_17net, highlight region 'SomMotA\_4 R', Scout function: mean



There should now be a new matrix object containing the time series data for the Schaefer parcel of interest. To export this time series for use in the burst detection tutorials:

- Right click on the matrix > File > Export to file > save as .mat file



For a simple transition to the Jupyter tutorials, save the matrix under the name 'sub-0002\_SomMotA\_4\_R.mat'

You should now be able to run any of the single channel burst detection tutorials including: Amplitude\_thresholding.ipynb, Cyclebycycle.ipynb, Microstates.ipynb, PAPTO.ipynb, HMM.ipynb, or eBOSC.ipynb.

## **5. A Note on Sensor-Level Multichannel Methods (CDL)**

For the sensor-level multichannel methods (i.e., CDL), extensive preprocessing in Brainstorm is not performed. For this method, source localization and time series extraction is not required. Therefore, for simplicity, the BIDS-formatted data is loaded directly into the Jupyter notebook and the preprocessing steps indicated in step 2 (e.g., filtering, artifact removal) are performed using MNE python prior to applying the method of interest. The preprocessing steps are matched to those used in the Brainstorm-processed data.

As such, CDL.ipynb can be run immediately following the download of the OMEGA sample data.